## Caso de Estudio 32

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- Akkocaoglu, M., Uysal, S., Tekdemir, I., Ake, K., et al. Implant design and interosseous stability of immediately placed implants: a human cadaver study. *Clinical Oral Implants Research*, 2005;16:202-209.
- Albrektsson, T., Brånemark, P.I., Hansson, H.A., Lindström J. Osseointegrated titanium implants. Requirements for ensuring a long-lasting, direct bone-to-implant anchorage in man. *Acta Orthop Scand*, 1981;52(2):155-70. Aparicio, C., Rangert, B., Sennerby, L. Immediate/early loading of dental implants: a report from the Sociedad Española de Implantes World Congress consensus meeting in Barcelona, Spain, 2002. *Clin Implant Dent Relat Res.*, 2003;5(1):57-60.
- Berglundh, T., Abrahamsson, I., Lang, N.P., Lindhe, J. De novo alveolar bone formation adjacent to endosseous implants. *Clin Oral Implants Res.* 2003;14(3):251-62. Botticelli, D., Persson, L.G., Lindhe, J., Berglundh, T. Bone tissue formation adjacent to implants placed in fresh extraction sockets: an experimental study in dogs. *Clin Oral Implants Res.* 2006;17(4):351-8.
- Botticelli, D., Renzi, A., Lindhe, J. and Berglundh, T., Implants in fresh extraction sockets: a prospective 5-year follow-up clinical study. *Clinical Oral Implants Research*, 2008;19:1226- 1232.
- Brånemark, P.I., Hansson, B.O., Adell, R., Breine, U., et al. Osseointegrated implants in the treatment of the edentulous jaw. Experience from a 10-year period. *Scandinavian Journal of Plastic and Reconstructive Surgery Supplementum*, 1997; 16:1-132.
- Brånemark, P.I., Osseointegration and its experimental background. *Journal of Prosthetic Dentistry*, 1983;50(3):399-410.
- Caneva, M., Salata, L.A., De Souza, S.S., Baffone, G., et al. Influence of implant positioning in extraction sockets on osseointegration: histomorphometric analyses in dogs. *Clin Oral Implants Res*, 2010a;21(1):43-9. Caneva, M., Salata, L.A., De Souza, S.S., Bressan, E., et al. Hard tissue formation adjacent to implants of various size and configuration immediately placed into extraction sockets: an experimental study in dogs. *Clinical Oral Implants Research*, 2010b;21: 885-890.
- Caneva, M., Botticelli, D., Viganò, P., Morellini, F., et al. Connective tissue grafts in conjunction with implants installed immediately into extraction sockets. An experimental study in dogs. *Clin Oral Implants Res.* 2013;24(1):50-6.
- Chen, S.T., Darby, I.B., Reynolds, E.C. A prospective clinical study of non-submerged immediate implants: clinical outcomes and esthetic results. *Clin Oral Implants Res.* 2007;18(5):552-62.
- Chen, S.T., Buser, D., Clinical and esthetic outcomes of implants placed in postextraction sites. *Int J Maxillofac Implants*, 2009;24:186-217.
- Cochran, D.L., Schenk, R.K., Lussi, A., Higginbottom, F.L., et al. Bone response to unloaded and loaded titanium implants with a sandblasted and acid-etched surface: a histometric study in the canine mandible. *J Biomed Mater Res.* 1998;40(1):1-11.

- Cochran, D.L., Buser, D., ten Bruggenkate, C.M., Weingart, D., et al. The use of reduced healing times on ITI implants with a sandblasted and acid-etched (SLA) surface: early results from clinical trials on ITI SLA implants. *Clin Oral Implants Res.* 2002;13(2):144-53.
- Cornelini, R., Rubini, C., Fioroni, M., Favero, G.A., et al. Transforming growth factor-beta 1 expression in the peri-implant soft tissues of healthy and failing dental implants. *J Periodonto.* 2003;74(4):446-50.
- Davies, J.E. In vitro modeling of the bone/implant interface. *Anat Rec.* 1996;245(2):426-45.
- Favero, G., Lang, N.P., De Santis, E., Gonzalez, G.B., et al. Ridge preservation at implants installed immediately after molar extraction. An experimental study in the dog. *Clin. Oral Impl. Res.* 2013;255-261.
- Ferrus, J., Cecchinato, D., Pjetursson, E.B., Lang, N.P., et al. Factors influencing ridge alterations following immediate implant placement into extraction sockets. *Clin Oral Implants Res.* 2010;21(1):22-9
- Gelb DA. Immediate implant surgery: three-year retrospective evaluation of 50 consecutive cases. *Int J Oral Maxillofac Implants.* 1993;8:388-399.
- Gomez-Roman, G., Schulte, W., d'Hoedt, B., Axman-Krcmar, D. The Frialit-2 implant, *Int J Oral Maxillofac Implants;* 1997;12 (3): 299-309.
- Horinchi, K., Uchida, H., Yamamoto, K., Sugimira, M. Immediate loading of Brånemark system implant following placement in edentulous patients: A clinical report. *Int J Oral Maxillofac Implants.* 15:824-830.
- Krump, J.L., Barnett, B.G., The immediate implant: a treatment alternative. *Int J Oral Maxillofac Implants.* 1991;6(1):19-23.
- Lang, N.P., Wetzel, A.C., Stich, H., Caffesse, R.G. Histologic probe penetration in healthy and inflamed peri-implant tissues. *Clin Oral Implants Res.* 1994;5(4):191-201.
- Lang, N.P., Berglundh, T., Heitz-Mayfield, L.J., Pjetursson, B.E., et al. Consensus statements and recommended clinical procedures regarding implant survival and complications. *Int J Oral Maxillofac Implants.* 2004;19:150-4.
- Nkenke, E., Hahn, M., Lell, M., Wiltfang, J., et al. Anatomic site evaluation of the zygomatic bone for dental implant placement. *Clinical Oral Implants Research,* 2003;14: 72-79.
- Quirynen, M., Van Assche, N., Botticelli, D., Berglundh, T. How does the timing of implant placement to extraction affect outcome? *Int J Maxiofac Implants,* 2007; 22: 202-203.
- Romanos, G., Toh, C.G., Siar, C.H., Swaminathan, D., et al. Peri-implant bone reactions to immediately loaded implants. An experimental study in monkeys. *J Periodontol,* 2001; 72(4): 506-511.
- Salvi, G.E., Lang, N.P. Diagnostic parameters for monitoring peri-implant conditions. *Int J Oral Maxillofac Implants.* 2004;19:116-27.
- Schnitman, P.A. Implant dentistry: where are we now? *J Am Dent Assoc.* 1993;124(4):38-47. Schulte, W., Heimke, G., Das Tubinger Sofort-implantat, *Quintessenz,* 1976;27(1).
- Schwartz-Arad D, Grossman Y, Chaushu G. The clinical effectiveness of implants placed immediately into fresh extraction sites of molar teeth. *J Periodontol* 2000; 71:839-844.
- Tarnow, D.P., Emtiaz, S., Classi, A. Immediate loading of threaded implants at stage 1 surgery in edentulous arches: ten consecutive case reports with 1- to 5-year data. *Int J Oral Maxillofac Implants.* 1997;12(3):319-24.
- Tomasi, C., Sanz, M., Cecchinato, D., Pjetursson, B., et al. Bone dimensional variations at implants placed in fresh extraction sockets: a multilevel multivariate analysis. *Clin Oral Implants Res.* 2010;21(1):30-6.
- Wilson, T.G. Jr, Schenk, R., Buser, D., Cochran, D. Implants placed in immediate extraction sites: a report of histologic and histometric analyses of human biopsies. *Int J Oral Maxillofac Implants.* 1998;13(3):333-41.

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