

Case study 61

Assisted Implant Rehabilitation in the Esthetic Zone Using TI-Base on Narrow Diameter Implants (NICE)

Dr. Stuardo Valenzuela Manfredi



Assisted implant rehabilitation in the esthetic zone using TI-Base on narrow diameter implants (NICE)

Author

Dr. Stuardo Valenzuela Manfredi received his DMD from Universidad Mayor, Santiago, Chile (2005).

In 2011, he completed his specialty in Oromaxillofacial Implantology at the Universidad de Chile (cum laude), and since 2012 he has been a member of the Oral Implantology Society, Chile. In 2015, Dr. Valenzuela received a Diploma in Digital Dentistry, Universidad de los Andes, Chile.

He is the Clinical Director of the Alpha-Bio Tec Chile Training Center, in which he is the leader speaker for introductory courses in surgery and prosthetics. Dr. Valenzuela manages VM Dental Studio, a private practice in Santiago, Chile, dedicated to implantology and aesthetic dentistry.

During his career, Dr. Valenzuela has participated in numerous rehabilitation and implant surgery courses and is a frequent guest speaker at industry seminars and conferences worldwide.



Abstract

Implants in the aesthetic zone is one of the main challenges in assisted implant rehabilitation. This clinical case present a tooth rehabilitation of young female patient with high esthetic expectations. The missing tooth was rehabilitated with Alpha-Bio Tec's NICE narrow diameter implants and screw-retained crowns using the new CAD/CAM TiBase.

Background

Implant-supported prostheses need to replicate both soft and hard tissues in order to achieve acceptable esthetics for patients.

After the teeth extraction, clinicians have several time points for implantation:

- 1. Immediate implantation in the fresh sockets following tooth extraction.
- 2. Early implantation 4 to 8 weeks following tooth extraction after soft tissue healing and after the mucosa is covering the socket entrance.
- 3. Delayed implantation 12 to 16 weeks following tooth extraction after substantial amount of new bone has formed.
- 4. Late implantation 6 months following tooth extraction after complete bone healing.

There are several indications for each treatment, depending on the initial clinical situation [2,3,4].

While implantation in fresh sockets with immediate implant loading has become more common in recent years, it is not always possible due to various factors such as [4,5]:

- Insufficient bone volume for dental implant stability
- Insufficient bone volume for the ideal 3D guided placement of our implants
- Oclussion
- Psychological profile

Additionally, the installation of dental implants in fresh sockets is associated with greater esthetic complications than implant installation in healed sockets. The most common of these complications is recession of the gingival margins (20% to 40%). This complication may be associated with the following factors [1]:

- Thin vestibular bone plate
- · Facial malposition of the implant
- Fine periodontal biotype

Early installation of implants after soft tissue healing is an optional treatment when there is no ideal situation for post-extraction sockets (presence of thick vestibular bone plate, absence of gingival recession, absence of infections). Early installation is performed after the soft tissue has healed, from 4 to 8 weeks after the extraction [6].

The objective of this treatment option is to have intact mucosa tissue during implant installation in order to perform a predictable Guided Bone Regeneration (GBR) in the vestibular area. This GBR is performed using a slow resorption collagen membrane and a graft material with a low resorption rate to compensate for changes in the ridge following tooth extraction [3].

The advantages of this technique are:

- Highly predictable esthetic results
- Lower risk of poor vestibular implant position
- Higher amounts of keratinized gingiva at the surgical site

Limitations / disadvantages of this technique are:

- Two surgeries, traumatic extraction without flap and implant installation with GBR
- A longer healing time is required when there are large periapical lesions

The right choice of the surgical treatment plan, and the restorative materials and abutment types, are the key to a successful outcome.

The increased demand for metal-free restorations in recent decades, has prompted the search for new materials with improved mechanical properties and esthetics. This, coupled with the introduction of CAD/CAM technology (computer-aided design/computer-aided manufacturing) in implantology, has begun to compete with the traditional way of making dental prostheses [9, 10].

Assisted implant rehabilitation in the esthetic zone using TI-Base on narrow diameter implants (NICE)

There are different types of abutments materials (metal or ceramic) and shapes (prefabricated or custom) available on the market. The material and the prosthesis type selected, have a high impact on the treatment plan, especially when working in the esthetic zone [8].

Carved prostheses using CAD/CAM processes have demonstrated a performance similar to those made by conventional methods, showing an acceptable marginal fit [7, 9, 10].

Using Ti-Base provides several benefits: a metal-free screw-retained crown in the esthetic zone, keeping the metal-to-metal prosthetic connection intact and avoiding the use of ceramics in the area where there is large concentration of stress. Short studies have demonstrated that the biomechanical and biological behavior are similar in both the conventional and the new CAD/CAM techniques [10].

Case overview

Female patient, 21 years of age, attends clinic consultation regarding tooth 22. The patient has no systemic diseases, is a nonsmoker, has no significant extraoral or intraoral findings.

Materials used

- NICE implant L11.5mm
- 1 Healing abutment, Ø 3.8 H5mm CHC
- 1 Esthetic simply straight abutment H 2mm CHC as a temporary abutment
- 1 CAD/CAM Titanium Base CHC
- 0.5cc of Alpha-Bio's GRAFT natural bovine bone (0.5-1.0mm)1cc
- 0.5cc Corticocancellous granules 0.5cc syringe
- AlphaBio's GRAFT Collagen membrane 20*30mm

Treatment plan

Following rigorous clinical examinations and complementary examinations (X-rays and CBCT), the decision is to install a narrow diameter 3.2mm implant (NICE, Alpha-Bio Tec.) plus horizontal guided bone regeneration of the defect with two layers of regeneration material (allograft and bovine xenograft), and a slow resorption collagen membrane Alpha-Bio's GRAFT, attached with titanium pins.

After 3 months, the implant connection provisionalization was performed, followed by the final rehabilitation with a fixed, screw-retained lithium disilicate prosthesis unit, and a Ti-Base (Alpha-Bio Tec).

Conclusion

This article describes a case of a single tooth implantation and restoration. As shown, the use of a narrow implant (NICE, Alpha-Bio Tec) and the CAD/CAM restoration technique enabled to receive a better aesthetic result for the patient.










Occlusal view of defect



3 Pre-surgery CBCT




4 Occlusal view of bone volume



5 Occlusal view of the implant installed. (NICE L11.5mm), CHC platform



6 GBR, in two layers of regenerative material, the first FDBA (AlphaBio's GRAFT) and a second layer of Bovine Xenograft (AlphaBio's GRAFT), plus a slow resorption collagen membrane (AlphaBio's GRAFT) attached with tacks and sutures



7 Post-surgery checkup CBCT



8 Connection of implant, occlusal view healing abutment



Provisionalization and connective graft (CTG) to improve gingival biotype



10 Provisional screwing of resin



11 Checkup 3 months following connection

Assisted implant rehabilitation in the esthetic zone using TI-Base on narrow diameter implants (NICE)



12 Occlusal view, emergency profile



Open tray printing, personalized transfer with acrylic resin



View of virtual working model and Ti-Base with 3Shape D750 Scanner



15 Crown Design



16 Ti-Base and lithium disilicate crown.



Lithium disilicate acid etching for 20 seconds with hydrofluoric acid



18 Lithium disilicate silanization for 60 seconds



19 Protection of prosthetic screw during cementing outside of the mouth



20 Cementing outside of the mouth with self-cured resin cement



21 Screw-retained crown cemented outside of the mouth




22 Checkup one month after installation



Bibliography

- 1. Araujo MG, Sukekava F, Wennstrom JL, Lindhe J. Ridge alterations following implant placement in fresh extraction sockets: An experimental study in the dog. J Clin Periodontol 2005;32:645-652.
- 2. Buser, Julia Wittneben Michael M. Bornstein, Linda Grutter, Vivianne Chappuis, and Urs C. Belser; Stability of Contour Augmentation and Esthetic Outcomes of Implant-Supported Single Crowns in the Esthetic Zone: 3-Year Results of a Prospective Study With Early Implant Placement Postextraction, J Periodonto, March 2011.
- 3. Buser D, Chen ST, Weber HP, Belser UC. Early implant placement following single-tooth extraction in the esthetic zone: Biologic rationale and surgical procedures. Int J Periodontics Restorative Dent 2008;28: 441-451.
- 4. Buser D, Chen ST. Implant placement in postextraction sites. In: Buser D, ed. 20 Years of Guided Bone Regeneration in Implant Dentistry, 2nd Edition. Chicago: Quintessence Publishing Co.; 2009:153-194.
- 5. Chen ST, Beagle J, Jensen SS, Chiapasco M, Darby I. Consensus statements and recommended clinical procedures regarding surgical techniques. Int J Oral Maxillofac Implants 2009;24 (Suppl):272-278.
- 6. Frank Higginbottom, Urs belser, Prostethic Management of Implants in the Estethics Zone, JOMI, Vol 19, 62-71, 2004.
- 7. Svanborg, Skjerven, Carlsson (2014), Marginal and Internal fit of cobalt chromium fixed dental prostheses generated from digital and conventional impressions. International Journal of Dentistry 2014.
- 8. William C. Martin, Robert A. Levine, Daniel Buser, Consensus Statements and Recommended Clinical Procedures Regarding Optimizing Esthetic Outcomes in Implant Dentistry, JOMI, Volume 29, Supplement, 2014.
- 9. Wismeijer, Urs Brägger Christopher Evans, ; Consensus Statements and Recommended Clinical Procedures Regarding Restorative Materials and Techniques for Implant Dentistry, JOMI, 2 Volume 28, Supplement, 2013.
- 10. Julia-Gabriela Wittneben, Robert F. Wright, Hans-Peter Weber, German O. Gallucci, ; A Systematic Review of the Clinical Performance of CAD/CAM Single-Tooth Restorations, The International Journal of prosthodontics, Vol 22, Number 5, 2009.
- 11. Theodoros Kapos, Christopher Evans; CAD/CAM Technology for Implants Abutments, Crowns, and superstructures, JOMI, Vol 29, 2014.







Smart Implantology Solutions

www. alpha-bio.net 💿 🕕 🔘





Alpha-Bio Tec Ltd.

7 Hatnufa St., P.O.B. 3936. Kiryat Arye, Petach Tikva 4951025, Israel T. +972.3.9291000 F. +972.3.9235055 sales@alpha-bio.net

MEDES LIMITED

5 Beaumont Gate, Shenley Hill Radlett, Herts WD7 7AR. England T./F. +44.192.3859810

International

T. +972.3.9291055 F. +972.3.9291010 export@alpha-bio.net